

Product datasheet

Anti-NOX3 antibody ab82708

★ ★ ★ ☆ ☆
[1 Abreviews](#)
[2 References](#)

Overview

Product name	Anti-NOX3 antibody
Description	Rabbit polyclonal to NOX3
Tested applications	Suitable for: IHC-Fr, WB
Species reactivity	Reacts with: Human Predicted to work with: Mouse
Immunogen	Synthetic peptide derived from an internal domain of Human NOX3

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: Whole serum
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab82708** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-Fr		1/50 - 1/500.
WB	★ ★ ★ ☆ ☆	1/200 - 1/2000. Predicted molecular weight: 65 kDa.

Target

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors