

Product datasheet

Natural Human D-Dimer protein ab98311

Overview

Product name	Natural Human D-Dimer protein
Protein length	Full length protein

Description

Nature	Native
Source	Native

Amino Acid Sequence

Species	Human
Molecular weight	200 kDa

Specifications

Our [Abpromise guarantee](#) covers the use of **ab98311** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Purity	> 90 % SDS-PAGE. Approximately 90% of the protein is D-Dimer as assessed by gradient polyacrylamide gel electrophoresis.
Form	Liquid
Additional notes	Plasma is donor tested and found negative for HIV I & II antibodies, Hepatitis B surface antigen, and Hepatitis C antibodies.

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Preservative: 0.09% Sodium Azide Constituents: 1% Human serum albumin, 0.15M Sodium chloride, 0.05M Monobasic dihydrogen sodium phosphate, pH 7.5
------------------------------	---

General Info

Relevance

A fibrin degradation fragment or product that is produced by the action of plasmin on fibrin in the clot dissolution process.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors