

Product datasheet

Natural human LDHA + LDHB protein **ab78637**

Overview

Product name	Natural human LDHA + LDHB protein
Protein length	Full length protein

Description

Nature	Native
Source	Native

Amino Acid Sequence

Species	Human
----------------	-------

Specifications

Our [Abpromise guarantee](#) covers the use of **ab78637** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Biological activity	Activity: 1060 Units/ml at 37°C. Specific Activity: 240 Units/mg Protein. One unit will catalyze the transformation of one micromole of L-lactate to pyruvate per minute at 37°C and pH 8.55.
----------------------------	---

Form	Liquid
-------------	--------

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C. Do Not Freeze. Preservative: None Constituents: 20mM Tris chloride, 1mM DTE, 1mM EDTA, pH 7.5 This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------------------------	--

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors