

Product datasheet

Normal Hamster Serum ab7483

1 Image

Overview

Product name	Normal Hamster Serum
Host species	Syrian hamster
Tested applications	Suitable for: Blocking, IHC-P, IHC-Fr
General notes	Total Protein Concentration : normal Strain: Golden Syrian Sex: Mixed sex Raised in: Syrian Hamster Purity: Whole antiserum

Properties

Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Purity	Whole antiserum

Applications

Our [Abpromise guarantee](#) covers the use of **ab7483** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Blocking		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

Normal Hamster Serum images



Other-Normal Hamster Serum(ab7483)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors