

Product datasheet

Anti-Phosphorylase B antibody ab31121

Overview

Product name	Anti-Phosphorylase B antibody
Description	Sheep polyclonal to Phosphorylase B
Host species	Sheep
Specificity	Recognizes Phosphorylase B. Exhibits high crossreactivity (95-135%) with Phosphorylase A by ELISA.
Tested applications	Suitable for: ELISA, WB
Species reactivity	Reacts with: Rabbit
Immunogen	Full length native protein: Phosphorylase B from rabbit muscle

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
Storage buffer	Preservative: None Constituents: Whole serum
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab31121** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/300 - 1/1000.
WB		1/500 - 1/1000. Predicted molecular weight: 121 kDa.

Target

Function	Phosphorylase b kinase catalyzes the phosphorylation of serine in certain substrates, including troponin I. The beta chain acts as a regulatory unit and modulates the activity of the holoenzyme in response to phosphorylation.
Pathway	Glycan biosynthesis; glycogen metabolism.
Involvement in disease	Glycogen storage disease 9B
Sequence similarities	Belongs to the phosphorylase b kinase regulatory chain family.
Post-translational modifications	Ser-701 is probably phosphorylated by PKA. Although the final Cys may be farnesylated, the terminal tripeptide is probably not removed, and the C-terminus is not methylated.
Cellular localization	Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors