

## Product datasheet

# Rabbit Anti-Cow IgA H&L (Alkaline Phosphatase) ab112627

### Overview

<b>Product name</b>	Rabbit Anti-Cow IgA H&L (Alkaline Phosphatase)
<b>Description</b>	Rabbit polyclonal Secondary Antibody to Cow IgA (AP)
<b>Host species</b>	Rabbit
<b>Target species</b>	Cow
<b>Tested applications</b>	<b>Suitable for:</b> WB, ELISA, IHC-P, ICC/IF
<b>Conjugation</b>	Alkaline Phosphatase

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	pH: 7.10 Preservative: 0.09% Sodium azide Constituents: 1.19% HEPES, 0.58% Sodium chloride, 0.01% Magnesium chloride, 0.001% Zinc chloride, 0.2% BSA
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	Antiserum was solid phase adsorbed to ensure class specificity. ab112627 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to alkaline phosphatase.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgA

### Applications

Our [Abpromise guarantee](#) covers the use of **ab112627** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000 - 1/10000. chemiluminescent detection.

Application	Abreviews	Notes
ELISA		1/1000 - 1/10000. as a primary.
IHC-P		1/100 - 1/1000.
ICC/IF		1/100 - 1/1000.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors