

Product datasheet

Rabbit Anti-Cow IgG H&L (HRP) ab6927

Overview

Product name	Rabbit Anti-Cow IgG H&L (HRP)
Target species	Cow
Tested applications	Suitable for: Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, WB, ICC/IF
Immunogen	Full length native bovine IgG (purified).
Conjugation	HRP

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.01% Gentamicin Sulfate Constituents: 10mg/ml BSA, 0.15M Sodium Chloride, 0.02M Potassium Phosphate. pH 7.2
Purity	IgG fraction
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using bovine IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
Conjugation notes	Horseradish Peroxidase (HRP)
Clonality	Polyclonal
General notes	This HRP conjugate antibody may have a small amount of sediments or cloudiness in the buffer.

Applications

Our [Abpromise guarantee](#) covers the use of **ab6927** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Dot blot		
ELISA		
IHC-P		

Application	Abreviews	Notes
IHC-Fr		
Immunomicroscopy		
WB		
ICC/IF		

Application notes

Dot: Use at an assay dependant dilution.

ELISA: 1/1000 - 1/5000.

IHC-P: Use at an assay dependant dilution.

IHC-Fr: Use at an assay dependant dilution.

IM: Use at an assay dependant dilution.

WB: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors