

Product datasheet

Rabbit Anti-Rat IgG H&L ab6703

7 References

Overview

Product name	Rabbit Anti-Rat IgG H&L
Target species	Rat
Tested applications	Suitable for: WB, ELISA, ICC/IF, IHC-P, IHC-Fr
Immunogen	Rat IgG whole molecule

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	0.02M K Phosphate, 0.15M NaCl, 0.01% Sodium Azide, pH7.2
Purity	IgG fraction
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rat IgG coupled to agarose beads.
Clonality	Polyclonal

Applications

Our [Abpromise guarantee](#) covers the use of **ab6703** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/2000 - 1/10000.
ELISA		1/20000 - 1/100000.
ICC/IF		1/1000 - 1/5000.
IHC-P		1/1000 - 1/5000.
IHC-Fr		1/1000 - 1/5000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors