

## Product datasheet

# Rabbit Anti-Sheep IgG H&L (Alkaline Phosphatase) ab6748

### 2 References

#### Overview

<b>Product name</b>	Rabbit Anti-Sheep IgG H&L (Alkaline Phosphatase)
<b>Host species</b>	Rabbit
<b>Target species</b>	Sheep
<b>Tested applications</b>	<b>Suitable for:</b> IHC-P, WB, ELISA, Immunomicroscopy, Dot blot, ICC/IF, IHC-Fr
<b>Immunogen</b>	Sheep IgG whole molecule
<b>Conjugation</b>	Alkaline Phosphatase

#### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C. Do Not Freeze.
<b>Storage buffer</b>	pH: 8 Preservative: 0.1% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.788% Tris HCl, 0.0095% Magnesium chloride, 0.014% Zinc chloride, 50% Glycerol
<b>Purity</b>	Affinity purified
<b>Purification notes</b>	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads
<b>Clonality</b>	Polyclonal
<b>General notes</b>	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons).

#### Applications

Our [Abpromise guarantee](#) covers the use of **ab6748** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		

Application	Abreviews	Notes
WB		
ELISA		
Immunomicroscopy		
Dot blot		
ICC/IF		
IHC-Fr		

**Application notes** suitable for immunoblotting (western or dot blot), ELISA and immunohistochemistry as well as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors