

Product datasheet

Recombinant Human CCL27 protein ab93953

1 Image

Overview

Product name Recombinant Human CCL27 protein
Protein length Full length protein

Description

Nature Recombinant
Source Escherichia coli

Amino Acid Sequence

Species Human
Sequence MFLPPSTAC CTQLYRKPLS DKLLRKVIQV
 ELQEADGDCH LQAFVLHLAQ RSICIHPQNP
 SLSQWFEHQE RKLHGTLPKL NFGMLRKMG

Specifications

Our [Abpromise guarantee](#) covers the use of **ab93953** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications Mass Spectrometry
 SDS-PAGE

Purity > 95 % SDS-PAGE.
 ab93953 was purified using conventional chromatography techniques.

Form Liquid

Additional notes Endotoxin Level: < 1.0 EU per 1 µg of protein (determined by LAL method)

Preparation and Storage


Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

Preservative: None
 Constituents: 10% Glycerol, 10mM Sodium citrate, pH 3.5

General Info

Function	Chemotactic factor that attracts skin-associated memory T-lymphocytes. May play a role in mediating homing of lymphocytes to cutaneous sites. Binds to CCR10.
Tissue specificity	Testis, thymus, placenta, ovary and skin.
Sequence similarities	Belongs to the intercrine beta (chemokine CC) family.
Cellular localization	Secreted.

Images



15% SDS-PAGE analysis of ab93953 (3µg)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors