

Product datasheet

Recombinant Human Dishevelled 3 protein ab114498

1 Image

Overview

Product name Recombinant Human Dishevelled 3 protein
Protein length Protein fragment

Description

Nature Recombinant
Source Wheat germ

Amino Acid Sequence

Accession [Q92997](#)
Species Human
Sequence MGETKIYHLDGQETPYLVKLPLPAERVTLADFKGVLQRPSYKFFFKSMD
 DDFGVVKEEISDDNAKLPCFNRRVSWLVSAEGSHDPAPFCADNPSELP
 *
Molecular weight 37 kDa including tags
Amino acids 1 to 100

Specifications

Our [Abpromise guarantee](#) covers the use of **ab114498** in the following tested applications.
 The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE
 Western blot
 ELISA
Form Liquid
Additional notes Protein concentration is above or equal to 0.05 mg/ml. This protein is best used within three months from the date of receipt.

Preparation and Storage

Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.


pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

General Info

Function	May play a role in the signal transduction pathway mediated by multiple Wnt genes.
Sequence similarities	Belongs to the DSH family. Contains 1 DEP domain. Contains 1 DIX domain. Contains 1 PDZ (DHR) domain.
Post-translational modifications	Ubiquitinated. Deubiquitinated by CYLD, which acts on 'Lys-63'-linked ubiquitin chains.
Cellular localization	Cytoplasm.

Images



ab114498 analysed by 12.5% SDS-PAGE and stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors