

Product datasheet

Recombinant Human Dishevelled 3 protein ab114498

1 Image

Overview

Product name	Recombinant Human Dishevelled 3 protein
Protein length	Protein fragment

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Accession	Q92997
Species	Human
Sequence	MGETKIIYHLDGQETPYLVKLPPLPAERVTLADFKGVLQRPSYKFFFKSMD DDFGVVKEEISDDNAKLPCFNRRVSVLWVSAEGSHDPDPAPFCADNPSELP *
Molecular weight	37 kDa including tags
Amino acids	1 to 100

Specifications

Our [Abpromise guarantee](#) covers the use of **ab114498** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Western blot ELISA
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml. This protein is best used within three months from the date of receipt.


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------------------------	---

General Info

Function	May play a role in the signal transduction pathway mediated by multiple Wnt genes.
Sequence similarities	Belongs to the DSH family. Contains 1 DEP domain. Contains 1 DIX domain. Contains 1 PDZ (DHR) domain.
Post-translational modifications	Ubiquitinated. Deubiquitinated by CYLD, which acts on 'Lys-63'-linked ubiquitin chains.
Cellular localization	Cytoplasm.

Recombinant Human Dishevelled 3 protein images



ab114498 analysed by 12.5% SDS-PAGE and stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors