

Product datasheet

Recombinant Human Eotaxin 3 protein ab87891

1 Image

Overview

Product name Recombinant Human Eotaxin 3 protein
Protein length Full length protein

Description

Nature Recombinant
Source Escherichia coli

Amino Acid Sequence

Species Human

Specifications

Our [Abpromise guarantee](#) covers the use of **ab87891** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE
Purity > 90 % SDS-PAGE.
 ab87891 was purified by conventional chromatography, after refolding of the isolated inclusion bodies in a renaturation buffer. Endotoxin level: < 1.0 EU per 1 microgram of protein (determined by LAL method).
Form Liquid

Preparation and Storage


Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
 Preservative: None
 Constituents: 10% Glycerol, PBS, pH 7.4

General Info

Function Chemotactic for eosinophils and basophils. Binds to CCR3.

Tissue specificity	Ubiquitously expressed at low levels in various tissues including heart and ovary.
Sequence similarities	Belongs to the intercrine beta (chemokine CC) family.
Cellular localization	Secreted.

Images



15% SDS-PAGE showing ab87891 (3 μ g).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors