

## Product datasheet

# Recombinant Human Flt3 / CD135 protein ab59307

### Overview

---

<b>Product name</b>	Recombinant Human Flt3 / CD135 protein
<b>Protein length</b>	Protein fragment

### Description

---

<b>Nature</b>	Recombinant
<b>Source</b>	Escherichia coli
<b>Amino Acid Sequence</b>	
<b>Species</b>	Human
<b>Sequence</b>	TQDCSFQHSPISDFAVKIRELSDYLLQDYPVTVASNLQDEELCGGLWRL VLA QRWMERLKTAVGSKMQGLLERVNTEIHFVTKCAFQPPPSCLRFVQ TNISRLQ ETSEQLVALKPWITRQNF SRCLELQCQPDSSTLPPWSPR PLEATAPTA
<b>Molecular weight</b>	18 kDa
<b>Amino acids</b>	27 to 181

### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab59307** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	SDS-PAGE HPLC
<b>Purity</b>	> 95 % SDS-PAGE. 0.2 µm sterile filtered. Endotoxin, tested. Suitable for cell culture.
<b>Form</b>	Lyophilised

### Preparation and Storage

---

<b>Stability and Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>Reconstitution</b>	Reconstitute in 5 to 50ul of distilled water (conc. 1 to 0.1mg/ml). This solution can then be diluted into other aqueous buffers and stored at 2-8 °C for up to 1 week. For extended storage, freeze in working aliquots.

## General Info

---

<b>Function</b>	Receptor for the FL cytokine. Has a tyrosine-protein kinase activity.
<b>Tissue specificity</b>	Bone marrow cells.
<b>Sequence similarities</b>	Belongs to the protein kinase superfamily. Tyr protein kinase family. CSF-1/PDGF receptor subfamily. Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Contains 1 protein kinase domain.
<b>Cellular localization</b>	Membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors