

Product datasheet

Recombinant Human Pleckstrin protein ab116861

1 Image

Overview

Product name Recombinant Human Pleckstrin protein

Protein length Protein fragment

Description

Nature Recombinant

Source Wheat germ

Amino Acid Sequence

Accession [P08567](#)

Species Human

Sequence PETIDLGALYLSMKDTEKGIKELNLEKDKKIFNHCFTGNCVIDWLVSQSVRNRQEGLMIASSLLNEGYLQPAGDMSKSAVDGTAENPFLDNPDAFYFPPDSGFFCEENS

Molecular weight 38 kDa including tags

Amino acids 121 to 230

Specifications

Our [Abpromise guarantee](#) covers the use of **ab116861** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications ELISA
SDS-PAGE
Western blot

Form Liquid

Additional notes Protein concentration is above or equal to 0.05 mg/ml.
Best use within three months from the date of receipt of this protein.

Preparation and Storage

Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

General Info

Function


Major protein kinase C substrate of platelets.

Sequence similarities

Contains 1 DEP domain.

Contains 2 PH domains.

Images



12.5% SDS-PAGE stained with Coomassie

Blue showing ab116861 at approximately

37.73 kDa.

SDS-PAGE - Pleckstrin protein (Human) (ab116861)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors