

Product datasheet

Recombinant Human RND3 protein ab97952

1 Image

Overview

Product name	Recombinant Human RND3 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Escherichia coli
Amino Acid Sequence	
Accession	P61587
Species	Human
Sequence	MGSSHHHHHSSGLVPRGSHMKERRASQKLSSKSIMDPNQNVKCKMVVG DSQCGKTALLHVFAKDCFPENYVPTVFENYASFEIDTQRIELSLWDTSG SPYYDNVRPLSYPDSDAVLICFDISRPETLDSVLKKWKGEIQEFCPNTKM LLVGCKSDLRTDVSTLVELSNHRQTPVSYDQGANNMAKQIGAATYIECSAL QSENSVRDIFHVATLACVNKTNKNVKNKSQRATKRISHMPSRPELSAVA TDLRKDKAKSC
Molecular weight	29 kDa including tags
Amino acids	1 to 241
Tags	His tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab97952** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Mass Spectrometry SDS-PAGE
Mass spectrometry	MALDI-TOF-TOF
Purity	> 90 % SDS-PAGE. ab97952 is purified using conventional chromatography techniques.
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


Preservative: None

Constituents: 40% Glycerol, 0.1M Sodium chloride, 20mM Tris HCl, 1mM EDTA, 1mM DTT, pH 8.0

General Info

Function	Binds GTP but lacks intrinsic GTPase activity and is resistant to Rho-specific GTPase-activating proteins.
Tissue specificity	Ubiquitous.
Sequence similarities	Belongs to the small GTPase superfamily. Rho family.
Cellular localization	Golgi apparatus membrane.

Recombinant Human RND3 protein images



15% SDS-PAGE showing ab97952 at approximately 29.2kDa (3µg).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors