

Product datasheet

Recombinant Human SF20 protein ab86915

1 Image

Overview

Product name	Recombinant Human SF20 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Escherichia coli

Amino Acid Sequence

Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MSEPTTVAFD VRPGGVVHSF SHNVGPGDKY TCMFTYASQG GTNEQWQMSL GTSEDHQHFT CTWRPQGKS LYFTQFKAE VRGAEIYAM AYSKAAFERE SDVPLKTEEF EVTKTAVAGR PGAFKAELSK LVVAKASRT EL

Specifications

Our [Abpromise guarantee](#) covers the use of **ab86915** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Purity	> 95 % SDS-PAGE. ab86915 was purified using conventional chromatography techniques, after refolding of the isolated inclusion bodies in a renaturation buffer.
Form	Liquid


Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Preservative: None Constituents: 20% Glycerol, 20mM Tris HCl, pH 8.0
------------------------------	--

General Info

Tissue specificity	Expressed in synovial tissue. Found in synovial fluid of patients with arthropaties.
Sequence similarities	Belongs to the UPF0556 family.
Cellular localization	Secreted. Endoplasmic reticulum-Golgi intermediate compartment.

Images



15% SDS-PAGE analysis of 3 ug ab86915

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors