

Product datasheet

Recombinant Human SF20 protein ab86915

1 Image

Overview

---

**Product name** Recombinant Human SF20 protein  
**Protein length** Full length protein

Description

---

**Nature** Recombinant  
**Source** Escherichia coli

Amino Acid Sequence

**Species** Human  
**Sequence** **MGSSHHHHH SSGLVPRGSH** MSEPTTVAFD  
 VRPGGVVHSF SHNVGPGDKY TCMFTYASQG  
 GTNEQWQMSL GTSEDHQHFT CTWRPQGKS  
 LYFTQFKAE VRGAEIYAM AYSKAAFERE  
 SDVPLKTEEF EVTKTAVAHR PGAFKAELSK  
 LVMKASRT EL

Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab86915** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** SDS-PAGE  
**Purity** > 95 % SDS-PAGE.  
 ab86915 was purified using conventional chromatography techniques, after refolding of the isolated inclusion bodies in a renaturation buffer.  
**Form** Liquid

Preparation and Storage

---

**Stability and Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.  
 Preservative: None  
 Constituents: 20% Glycerol, 20mM Tris HCl, pH 8.0

## General Info


---

<b>Tissue specificity</b>	Expressed in synovial tissue. Found in synovial fluid of patients with arthropaties.
<b>Sequence similarities</b>	Belongs to the UPF0556 family.
<b>Cellular localization</b>	Secreted. Endoplasmic reticulum-Golgi intermediate compartment.

---

## Images

---



15% SDS-PAGE analysis of 3 ug ab86915

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors