

Product datasheet

Recombinant Human ZNF207 protein ab159878

1 Image

Overview

<b>Product name</b>	Recombinant Human ZNF207 protein
<b>Protein length</b>	Full length protein

Description

<b>Nature</b>	Recombinant
<b>Source</b>	Wheat germ

Amino Acid Sequence

<b>Species</b>	Human
<b>Sequence</b>	MGRK K K K K Q L K P W C W Y C N R D F D D E K I L I Q H Q K A K H F K C H I C H K K L Y T G P G L A I H C M Q V H K E T I D A V P N A I P G R T D I E L E Y G M E G I P E K D M D E R R R L L E Q K T Q E S Q K K K Q Q D D S D E Y D D D D S A A S T S F Q P Q P V Q P Q Q G Y I P P M A Q P G L P P V P G A P G M P P G I P L M P G V P P L M P G M P P V M P G M P P G L H H Q R K Y T Q S F C G E N I M M P M G G M M P P G P I P L M P G M P P G M P P V P R P G I P P M T Q A Q A V S A P G I L N R P P A P T A T V P A P Q P P V T K P L F P S A G Q A Q A A V Q G P V G T D F K P L N S T P A T T T E P P K P T F P A Y T Q S T A S T T S T T N S T A A K P A A S I T S K P A T L T T T S A T S K L I H P D E D I S L E E R R A Q L P K Y Q R N L P R P G Q A P I G N P P V G P I G G M M P P Q P G I P Q Q Q G M R P P M P P H G Q Y G G H H Q G M P G Y L P G A M P P Y G Q G P P M V P P Y Q G G P P R P P M G M R P P V M S Q G G R Y

<b>Amino acids</b>	1 to 463
<b>Tags</b>	proprietary tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab159878** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	Western blot ELISA
<b>Form</b>	Liquid
<b>Additional notes</b>	Protein concentration is above or equal to 0.05 mg/ml.

Preparation and Storage

## Preparation and Storage

### Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

## General Info

### Tissue specificity

Ubiquitous.

### Sequence similarities

Contains 2 C2H2-type zinc fingers.

### Cellular localization

Nucleus.

## Images



ab159878 on a 12.5% SDS-PAGE stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors