


Product datasheet

Anti-Resistin antibody ab93069

1 Image

Overview

<b>Product name</b>	Anti-Resistin antibody
<b>Description</b>	Chicken polyclonal to Resistin
<b>Host species</b>	Chicken
<b>Tested applications</b>	<b>Suitable for:</b> WB, Indirect ELISA
<b>Species reactivity</b>	<b>Reacts with:</b> Recombinant fragment <b>Predicted to work with:</b> Mouse, Rat  <b>Does not react with:</b> Human
<b>Immunogen</b>	Synthetic peptide: VSSRGRLASCPEGT (Mouse)  <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
<b>Storage buffer</b>	Preservative: 0.02% Sodium Azide Constituents: PBS
<b>Purity</b>	Protein L purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgY

Applications

Our [Abpromise guarantee](#) covers the use of **ab93069** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		
Indirect ELISA		

Application	Abreviews	Notes
-------------	-----------	-------

#### Application notes

I-ELISA: 1/2000 - 1/5000.

WB: 1/2000 - 1/5000. Predicted molecular weight: 11 kDa.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

#### Target

##### Function

Hormone that seems to suppress insulin ability to stimulate glucose uptake into adipose cells. Potentially links obesity to diabetes.

##### Tissue specificity

Expressed only in fatty tissues.


##### Sequence similarities

Belongs to the resistin/FIZZ family.

##### Cellular localization

Secreted.

#### Images



Anti-Resistin antibody (ab93069) at 1/2000

dilution + recombinant Mouse Resistin at 0.02

µg

**Predicted band size:** 11 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors