


Product datasheet

Anti-Resistin antibody ab93069

1 Image

Overview

Product name	Anti-Resistin antibody
Description	Chicken polyclonal to Resistin
Host species	Chicken
Tested applications	Suitable for: WB, Indirect ELISA
Species reactivity	Reacts with: Recombinant fragment Predicted to work with: Mouse, Rat  Does not react with: Human
Immunogen	Synthetic peptide: VSSRGRLASCPEGT (Mouse) Run BLAST with Run BLAST with

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: 0.02% Sodium Azide Constituents: PBS
Purity	Protein L purified
Clonality	Polyclonal
Isotype	IgY

Applications

Our [Abpromise guarantee](#) covers the use of **ab93069** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		
Indirect ELISA		

Application	Abreviews	Notes
-------------	-----------	-------

Application notes

I-ELISA: 1/2000 - 1/5000.
 WB: 1/2000 - 1/5000. Predicted molecular weight: 11 kDa.

Not yet tested in other applications.
 Optimal dilutions/concentrations should be determined by the end user.

Target

Function

Hormone that seems to suppress insulin ability to stimulate glucose uptake into adipose cells. Potentially links obesity to diabetes.

Tissue specificity

Expressed only in fatty tissues.


Sequence similarities

Belongs to the resistin/FIZZ family.

Cellular localization

Secreted.

Images



Anti-Resistin antibody (ab93069) at 1/2000 dilution + recombinant Mouse Resistin at 0.02 µg

Predicted band size: 11 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors