

## Product datasheet

# Anti-SESN2 antibody ab115984

1 Image

### Overview

|                            |   |
|----------------------------|---|
| <b>Product name</b>        | Anti-SESN2 antibody   |
| <b>Description</b>         | Rabbit polyclonal to SESN2  |
| <b>Tested applications</b> | <b>Suitable for:</b> WB   |
| <b>Species reactivity</b>  | <b>Reacts with:</b> Human<br><b>Predicted to work with:</b> Mouse, Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog, Pig   |
| <b>Immunogen</b>           | Synthetic peptide corresponding to a region within C terminal amino acids 268-317 (LQESLLRDEG TSQEEMESRF ELEKSESLLV TPSADILEPS PHPDMLCFVE) of Human SESN2 (NP_113647).<br><a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a> |
| <b>Positive control</b>    | Human fetal heart lysate.   |

### Properties

|                             |  |
|-----------------------------|--|
| <b>Form</b>                 | Liquid   |
| <b>Storage instructions</b> | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| <b>Storage buffer</b>       | Constituents: 98% PBS, 2% Sucrose  |
| <b>Purity</b>               | Immunogen affinity purified  |
| <b>Clonality</b>            | Polyclonal   |
| <b>Isotype</b>              | IgG  |

### Applications

Our [Abpromise guarantee](#) covers the use of **ab115984** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| Application | Abreviews | Notes  |
|-------------|-----------|--|
| WB          |           | Use a concentration of 1 µg/ml. Predicted molecular weight: 54 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |

## Target

---

## Images

---



Anti-SESN2 antibody (ab115984) at 1 µg/ml +  
Human fetal heart lysate at 10 µg

**Predicted band size : 54 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors