

Product datasheet

SIVA overexpression 293T lysate (whole cell) ab94200

[2 Images](#)

Overview

Product name	SIVA overexpression 293T lysate (whole cell)
Description	SIVA overexpression lysate
General notes	ab94200 is a 293T cell transfected lysate in which Human SIVA has been transiently over-expressed using a pCMV-SIVA plasmid. The lysate is provided in 1X Sample Buffer. Note: For more detailed how the transfected lysate was prepared view preparation notes
Tested applications	Suitable for: WB

Properties

Form	Liquid
Storage instructions	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	pH: 6.80 Constituent: 100% 1x Sample Buffer
Purity	Whole Cell Lysate
Background	Function: Induces CD27-mediated apoptosis. Inhibits BCL2L1 isoform Bcl-x(L) anti-apoptotic activity. Inhibits activation of NF-kappa-B and promotes T-cell receptor-mediated apoptosis. Tissue specificity: Ubiquitous. Mostly expressed in thymus, testis, ovary, prostate, small intestine and spleen and less in colon. PTM: Phosphorylated by ABL2/ARG in response to oxidative stress.


Applications

Our [Abpromise guarantee](#) covers the use of **ab94200** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent dilution.

Images



ab94200 at 15µg/lane on an SDS-PAGE gel.



All lanes : Anti-SIVA antibody (ab67620) at
1/500 dilution

Lane 1 : SIVA overexpression 293T lysate
(whole cell) (ab94200)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

Secondary

Goat Anti-mouse IgG (H and L) HRP conjug at
1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors