

Product datasheet

Sheep Anti-Horse IgA H&L ab112881

Overview

| | |
|----------------------------|--|
| Product name | Sheep Anti-Horse IgA H&L |
| Description | Sheep polyclonal Secondary Antibody to Horse IgA |
| Host species | Sheep |
| Target species | Horse |
| Tested applications | Suitable for: WB, ELISA, ICC/IF, IHC-P |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| Storage buffer | pH: 7.20 Preservative: 0.1% Sodium azide Constituent: 99% PBS |
| Purity | Immunogen affinity purified |
| Purification notes | Antiserum was solid phase adsorbed to ensure subclass specificity. ab112881 was isolated by affinity chromatography using antigen coupled to agarose beads. |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications

Our [Abpromise guarantee](#) covers the use of **ab112881** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-------------|-----------|--|
| WB | | 1/1000 - 1/30000. Colorimetric: Use 1/1000-1/10000 dilution; Chemiluminescent: Use 1/1000-1/30000 dilution |
| ELISA | | 1/100 - 1/30000. Coating: Use 1/100-1/500 dilution; Primary: Use 1/1000-1/30000 dilution. |
| ICC/IF | | 1/200 - 1/2000. |

| Application | Abreviews | Notes |
|-------------|-----------|-----------------|
| IHC-P | | 1/200 - 1/2000. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors