

## Product datasheet

# Sheep Anti-Rabbit IgG H&L (HRP) ab97095

### 1 References

#### Overview

<b>Product name</b>	Sheep Anti-Rabbit IgG H&L (HRP)
<b>Host species</b>	Sheep
<b>Target species</b>	Rabbit
<b>Specificity</b>	By immunoelectrophoresis and ELISA this antibody reacts specifically with Rabbit IgG and with light chains common to other Rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins.
<b>Tested applications</b>	<b>Suitable for:</b> ICC, IHC-P, ELISA, WB
<b>Conjugation</b>	HRP

#### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	Preservative: 0.1% Proclin Constituents: 0.2% BSA, PBS
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Horse Radish Peroxidase (HRP).
<b>Conjugation notes</b>	Molar enzyme: antibody protein ratio is 1:1
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

#### Applications

Our [Abpromise guarantee](#) covers the use of **ab97095** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC		Use at an assay dependent dilution.
IHC-P		1/200 - 1/5000.

Application	Abreviews	Notes
ELISA		1/10000 - 1/100000. (Primary ELISA)
WB		1/5000 - 1/50000. (Colorimetric: 1/5000-1/30000; Chemiluminescent: 1/10000-1/50000)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors