

Product datasheet

Anti-Urotensin II antibody (Biotin) ab48319

Overview

<b>Product name</b>	Anti-Urotensin II antibody (Biotin)
<b>Description</b>	Rabbit polyclonal to Urotensin II (Biotin)
<b>Host species</b>	Rabbit
<b>Conjugation</b>	Biotin
<b>Specificity</b>	This antibody is specific for Urotensin II.
<b>Tested applications</b>	<b>Suitable for:</b> RIA, ELISA, IP
<b>Species reactivity</b>	<b>Reacts with:</b> Human
<b>Immunogen</b>	Synthetic peptide: AGTADCFWKYCV conjugated to KLH via carboxyl group

 [Run BLAST with](#)

 [Run BLAST with](#)

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 50% Glycerol, PBS, pH 7.5
<b>Purity</b>	Protein G purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab48319** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
RIA		Use at an assay dependent dilution.

Application	Abreviews	Notes
ELISA		Use at an assay dependent dilution.
IP		Use at an assay dependent dilution.

## Target

<b>Function</b>	Highly potent vasoconstrictor.
<b>Tissue specificity</b>	Brain specific.
<b>Sequence similarities</b>	Belongs to the urotensin-2 family.
<b>Cellular localization</b>	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors