

## Product datasheet

# Anti-VSV-G tag antibody ab3861

### 3 References

#### Overview

<b>Product name</b>	Anti-VSV-G tag antibody
<b>Description</b>	Goat polyclonal to VSV-G tag
<b>Host species</b>	Goat
<b>Tested applications</b>	<b>Suitable for:</b> IP, ELISA, ICC, WB
<b>Immunogen</b>	Synthetic peptide: YTDIEMNRLGK conjugated to KLH.

 [Run BLAST with](#)

 [Run BLAST with](#)

#### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>Storage buffer</b>	PBS with 0.1% sodium azide pH 7.2
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

#### Applications

Our [Abpromise guarantee](#) covers the use of **ab3861** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IP		Use at an assay dependent dilution.
ELISA		1/100 - 1/500. For coating. Primary: 1/1000 – 1/30000.
ICC		1/200 - 1/2000.

Application	Abreviews	Notes
WB		1/1000 - 1/30000.

---

## Target

**Relevance** Vesicular stomatitis virus (VSV), an enveloped RNA virus from the Rhabdoviridae family, is released from the plasma membrane of host cells by a process called budding. The glycoprotein (VSV-g) contains a domain in its extracellular membrane proximal stem that appears to be needed for efficient VSV budding.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

---

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors