

Product datasheet

Anti-Vibrio parahaemolyticus O8 antibody ab79632

Overview

Product name	Anti-Vibrio parahaemolyticus O8 antibody
Description	Rabbit polyclonal to Vibrio parahaemolyticus O8
Host species	Rabbit
Specificity	Specific for Vibrio parahaemolyticus somatic O8 antigen
Tested applications	Suitable for: Agglutination
Species reactivity	Vibrio parahaemolyticus O8
Immunogen	Tissue/ cell preparation of Vibrio parahaemolyticus O8 (heat inactivated)
Positive control	Antigenic suspension.
General notes	Reagents should be allowed to stand at 15°C - 25°C for at least 30 minutes prior to use.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.1% Sodium Azide Constituents: Whole serum
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab79632** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Agglutination		Use at an assay dependent concentration. See attached protocol.

Target

Relevance

Vibrio parahaemolyticus is a curved, rod-shaped, Gram-negative bacterium found in brackish saltwater, which, when ingested, causes gastrointestinal illness in humans. *V. parahaemolyticus* is oxidase positive, facultatively aerobic, and does not form spores. Like other members of the genus *Vibrio*, this species is motile, with a single, polar flagellum.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors