


Product datasheet

Anti-ZCCHC7 antibody ab104503

2 Images

Overview

<b>Product name</b>	Anti-ZCCHC7 antibody
<b>Description</b>	Rabbit polyclonal to ZCCHC7
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> WB, IHC-P
<b>Species reactivity</b>	<b>Reacts with:</b> Human <b>Predicted to work with:</b> Mouse, Rat 
<b>Immunogen</b>	Recombinant fragment, corresponding to amino acids 83-332 of Human ZCCHC7 (BC036940).
<b>Positive control</b>	293T cell lysate; Human Fetal Brain lysate; Human small intestine

Properties

<b>Form</b>	Lyophilised:Reconstitute in 200ul sterile H2O.
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
<b>Storage buffer</b>	Preservative: 0.02% Sodium Azide Constituents: PBS, pH 7.2
<b>Purity</b>	Protein A purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab104503** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/1000. Predicted molecular weight: 63 kDa.
IHC-P		1/100 - 1/500.

## Target

<b>Sequence similarities</b>	Contains 4 CCHC-type zinc fingers.
<b>Cellular localization</b>	Nucleus, nucleolus.

## Images




**All lanes :** Anti-ZCCHC7 antibody  
(ab104503) at 1/500 dilution

**Lane 1 :** 293T cell lysate

**Lane 2 :** Human fetal brain lysate

**Predicted band size:** 63 kDa



ab104503, at 1/100 dilution, staining  
ZCCHC7 in formalin-fixed, paraffin-embedded  
Human small intestine by  
Immunohistochemistry.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors