

Product datasheet

Anti-BMP8b antibody ab110048

1 Image

Overview

Product name	Anti-BMP8b antibody
Description	Rabbit polyclonal to BMP8b
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Synthetic peptide derived from the internal sequence of Human BMP8
Positive control	Extracts from 293 and A549 cells

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
Storage buffer	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: PBS, 0.88% Sodium chloride, 50% Glycerol
Purity	Immunogen affinity purified
Purification notes	ab110048 is affinity-purified from rabbit antiserum by affinity-chromatography using an epitope-specific immunogen.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab110048** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/1000. Predicted molecular weight: 45 kDa.

Target

Function	Induces cartilage and bone formation. May be the osteoinductive factor responsible for the phenomenon of epithelial osteogenesis. Plays a role in calcium regulation and bone homeostasis.
Sequence similarities	Belongs to the TGF-beta family.
Cellular localization	Secreted.

Images



All lanes : Anti-BMP8b antibody (ab110048)
at 1/500 dilution

Lane 1 : extracts from A549 cells

Lane 2 : extracts from 293 cells

Lane 3 : extracts from A549 cells with
blocking peptide at 10 µg

Lysates/proteins at 30 µg per lane.

Predicted band size: 45 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors