


Product datasheet

Anti-ZBT24 antibody ab111404

1 Image

Overview

Product name	Anti-ZBT24 antibody
Description	Rabbit polyclonal to ZBT24
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Mouse Predicted to work with: Human 
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 144-174 of Human ZBT24 (NP_055612.1 , O43167)
Positive control	Mouse cerebellum lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
Storage buffer	Preservative: 0.09% Sodium azide Constituent: PBS
Purity	Immunogen affinity purified
Purification notes	Purified through a protein A column, followed by peptide affinity purification
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab111404** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/100 - 1/500. Predicted molecular weight: 78 kDa.

Target

Function	May be involved in BMP2-induced transcription.
Sequence similarities	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 1 A.T hook DNA-binding domain. Contains 1 BTB (POZ) domain. Contains 8 C2H2-type zinc fingers.
Cellular localization	Nucleus.

Images



Anti-ZBT24 antibody (ab111404) at 1/100 dilution + Mouse cerebellum tissue lysate at 35 µg

Predicted band size : 78 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors