

Product datasheet

Anti-FeSODB1 antibody ab115977

Overview

Product name	Anti-FeSODB1 antibody
Description	Chicken polyclonal to FeSODB1
Host species	Chicken
Tested applications	Suitable for: WB, ELISA
Species reactivity	Reacts with: Leishmania
Immunogen	Synthetic peptide (Leishmania) Purified from FeSODB1 protein, covalently modified to carrier protein and emulsified in adjuvant.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: 0.02% Sodium azide
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgY

Applications

Our [Abpromise guarantee](#) covers the use of **ab115977** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500. Predicted molecular weight: 22 kDa.
ELISA		1/500.

Target

Relevance	Destroys radicals which are normally produced within the cells and which are toxic to biological
------------------	--

systems.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors