

Product datasheet

Anti-WDR87 antibody ab115981

1 Image

Overview

| | |
|----------------------------|--|
| Product name | Anti-WDR87 antibody |
| Description | Rabbit polyclonal to WDR87 |
| Host species | Rabbit |
| Tested applications | Suitable for: WB |
| Species reactivity | Reacts with: Human |
| Immunogen | Synthetic peptide corresponding to a region within C terminal amino acids 810-859 (PQRELEWDRS QEFFFWHSRV RAISNTEYPK NKEEDEHFLE MRLSKDVTYS) of Human WDR87 (BAC87627). Run BLAST with ExPASy Run BLAST with NCBI |
| Positive control | A549 cell lysate. |

Properties

| | |
|-----------------------------|--|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| Storage buffer | Constituents: 98% PBS, 2% Sucrose |
| Purity | Immunogen affinity purified |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications

Our [Abpromise guarantee](#) covers the use of **ab115981** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| Application | Abreviews | Notes |
|-------------|-----------|---|
| WB | | Use a concentration of 1 µg/ml. Predicted molecular weight: 150 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |

Target

Sequence similarities

Contains 7 WD repeats.

Images



Anti-WDR87 antibody (ab115981) at 1 µg/ml
+ A549 cell lysate at 10 µg

Predicted band size: 150 kDa

Gel concentration: 6-18%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors