

Product datasheet

Anti-Meis homeobox 3 antibody ab116003

1 Image

Overview

Product name	Anti-Meis homeobox 3 antibody
Description	Rabbit polyclonal to Meis homeobox 3
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Rat Predicted to work with: Mouse, Rabbit, Horse, Guinea pig, Cow, Cat, Dog, Human
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 301-350 (LAQDTGLTIL QVNNWFINAR RRVQPMIDQ SNRTGQGASF NPEGQSMAGY) of Rat Meis homeobox 3(NP_001101942). Run BLAST with ExPASy Run BLAST with NCBI
Positive control	Rat lung lysate.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Constituents: 98% PBS, 2% Sucrose
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab116003** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 45 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

Target

Sequence similarities	Belongs to the TALE/MEIS homeobox family. Contains 1 homeobox DNA-binding domain.
Cellular localization	Nucleus.

Images



Anti-Meis homeobox 3 antibody (ab116003)
at 1 µg/ml + Rat lung lysate at 10 µg

Predicted band size: 45 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors