

Product datasheet

Anti-COMT antibody ab118908

2 Images

Overview

Product name	Anti-COMT antibody
Description	Rabbit polyclonal to COMT
Tested applications	Suitable for: IHC-P
Species reactivity	Reacts with: Human
Immunogen	Synthetic peptide conjugated to KLH, corresponding to the N-terminal domain of Human COMT.
Positive control	Human Adrenal and Human Liver tissues

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
Storage buffer	Preservative: 0.1% Sodium azide Constituent: 99% PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab118908** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use a concentration of 5 µg/ml.

Target

Function	Catalyzes the O-methylation, and thereby the inactivation, of catecholamine neurotransmitters
-----------------	---

and catechol hormones. Also shortens the biological half-lives of certain neuroactive drugs, like L-DOPA, alpha-methyl DOPA and isoproterenol.

Tissue specificity

Brain, liver, placenta, lymphocytes and erythrocytes.

Sequence similarities

Belongs to the mammalian catechol-O-methyltransferase family.


Post-translational modifications

The N-terminus is blocked.

Cellular localization


Cytoplasm and Cell membrane.

Images



Human Adrenal tissue, Formalin-Fixed, Paraffin-Embedded (FFPE) stained using ab118908 at a concentration of 5 µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-COMT antibody (ab118908)



Human Liver tissue Formalin-Fixed, Paraffin-Embedded (FFPE) stained using ab118908 at a concentration of 5 µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-COMT antibody (ab118908)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors