

Product datasheet

Anti-NKAIN2 antibody ab126547

1 Image

Overview

Product name	Anti-NKAIN2 antibody
Description	Rabbit polyclonal to NKAIN2
Host species	Rabbit
Tested applications	Suitable for: IHC-P
Species reactivity	Reacts with: Human
Immunogen	antigen sequence: GCTVTSVTPA PDWAPEDHRY ITVSGCLLEY QYIEVAH, corresponding to amino acids 114-150 of Human NKAIN2. Run BLAST with ExPASy Run BLAST with NCBI
Positive control	Human placenta tissue.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	pH: 7,20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab126547** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/20 - 1/50. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

Target

Tissue specificity	Expressed in fetal brain. Weakly expressed in adult brain and thymus. Not expressed in any other normal tissue examined.
Involvement in disease	Note=A chromosomal aberration involving NKAIN2 is a cause of lymphoma. Deletion del(6)(q13q21) within NKAIN2 and involving SUSP1 generates the SUSP1-NKAIN2 product.
Sequence similarities	Belongs to the NKAIN family.
Cellular localization	Cell membrane.

Images



ab126547 at 1/20 dilution staining NKAIN2 in paraffin embedded Human placenta tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NKAIN2 antibody (ab126547)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors