

Product datasheet

Anti-Bacillus anthracis protective antigen antibody [BAP0104] ab1991

1 References

Overview

Product name	Anti-Bacillus anthracis protective antigen antibody [BAP0104]
Description	Mouse monoclonal [BAP0104] to Bacillus anthracis protective antigen
Host species	Mouse
Specificity	This antibody is not cross-reacting with lethal factor antigen of B. anthracis, Y.Pestis, F. Tularensis, and Toxoplasma gondii.
Tested applications	Suitable for: ELISA
Immunogen	Highly purified protective antigen of Bacillus anthracis.
General notes	Concentration varies from lot to lot and can be provided on request.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C. Do Not Freeze.
Storage buffer	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.4
Purity	Protein G purified
Purification notes	Purity is tested by electrophoresis.
Clonality	Monoclonal
Clone number	BAP0104
Myeloma	Sp2/0
Isotype	IgG2b

Applications

Our [Abpromise guarantee](#) covers the use of **ab1991** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent dilution.

Target

Relevance	Anthrax toxins are composed of three distinct proteins, a protective antigen (PA), a lethal factor (LF) and an edema factor (EF). None of these proteins are toxic by themselves and several studies indicate that the anthrax toxin has the familiar A-B enzymatic binding structure, with PA acting as the binding domain and EF and/or LF acting as the active fragments.
Cellular localization	Secreted

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors