


Product datasheet

Anti-DMRT8.2 antibody ab36912

Overview

| | |
|---------------------------|---|
| Product name | Anti-DMRT8.2 antibody |
| Description | Rabbit polyclonal to DMRT8.2 |
| Host species | Rabbit |
| Species reactivity | Predicted to work with: Mouse  |
| Immunogen | Synthetic peptide conjugated to KLH derived from within residues 1 - 100 of Mouse DMRT8.2. Read Abcam's proprietary immunogen policy (Peptide available as ab36911 .) |

Properties

| | |
|-----------------------------|--|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| Storage buffer | Preservative: 0.02% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4 |
| Purity | Immunogen affinity purified |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications

| | |
|--------------------------|---|
| Application notes | This antibody gave a positive result in ELISA against the immunizing peptide (ab36911). Not yet tested in other applications. |
|--------------------------|---|

Target

| | |
|------------------------------|--|
| Relevance | The vertebrate DMRT gene family encodes putative transcription factors related to the sexual regulators Doublesex (Drosophila) and MAB-3 (Worm). DMRT8 is related to DMRT7 and located on the X chromosome in placental mammals. DMRT8.2 expression is restricted to male and female gonads and may be involved in mouse sexual development. |
| Cellular localization | Nuclear |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors