

## Product datasheet

# Adipose (human) normal tissue slides ab4274

### Overview

<b>Product name</b>	Adipose (human) normal tissue slides
<b>General notes</b>	<p>After excision, tissue was fixed in paraformaldehyde and processed to paraffin wax. Tissue sections were cut to a thickness of 5µm and mounted on positively charged glass slides. One tissue from each lot was stained by haematoxylin and eosin (H+E) and then microscopically examined to ensure quality. All tissue is examined / diagnosed by a licensed pathologist and is ethically obtained. Sex and age of the donor is lot dependent. Please enquire.</p> <p>This product is an excellent control for use in immunohistochemistry (IHC) and in-situ hybridization.</p>
<b>Tested applications</b>	<b>Suitable for:</b> IHC-P, In situ hybridization

### Properties

<b>Storage instructions</b>	Store at +4°C.
<b>Slide notes</b>	This product is an excellent control for use in immunohistochemistry (IHC) and in-situ hybridization.

### Applications

Our [Abpromise guarantee](#) covers the use of **ab4274** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use at an assay dependent concentration.
In situ hybridization		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors