

Product datasheet

Artery (human) normal tissue slides ab4283

Overview

Product name	Artery (human) normal tissue slides
General notes	<p>After excision, tissue was fixed in paraformaldehyde and processed to paraffin wax. Tissue sections were cut to a thickness of 5µm and mounted on positively charged glass slides. One tissue from each lot was stained by haematoxylin and eosin (H+E) and then microscopically examined to ensure quality. All tissue is examined / diagnosed by a licensed pathologist and is ethically obtained. Sex and age of the donor is lot dependent. Please enquire.</p> <p>This product is an excellent control for use in immunohistochemistry (IHC) and in-situ hybridization.</p>
Tested applications	Suitable for: IHC-P, In situ hybridization

Properties

Storage instructions	Store at +4°C.
Slide notes	This product is an excellent control for use in immunohistochemistry (IHC) and in-situ hybridization.

Applications

Our [Abpromise guarantee](#) covers the use of **ab4283** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use at an assay dependent concentration.
In situ hybridization		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors