

Product datasheet

ADAM29 peptide ab46842

Overview

Product name ADAM29 peptide

Description

Nature Synthetic

Specifications

Our [Abpromise guarantee](#) covers the use of **ab46842** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications ELISA
Western blot
Blocking

Purity > 95 % SDS-PAGE.

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Preservative: None
Constituents: 0.001% Tween 20, PBS, 30mM HEPES, 2mM EDTA, 150mM Sodium chloride, pH 6.75

General Info

Function May be involved in spermatogenesis and fertilization. Seems to be a non catalytic metalloprotease-like protein.

Tissue specificity Expressed specifically in testes.

Sequence similarities Contains 1 disintegrin domain.
Contains 1 EGF-like domain.

Contains 1 peptidase M12B domain.

Cellular localization

Membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors