

Product datasheet

Anti-LAMP1 antibody [B-T47] (Phycoerythrin) ab47210

Overview

Product name	Anti-LAMP1 antibody [B-T47] (Phycoerythrin)
Description	Mouse monoclonal [B-T47] to LAMP1 (Phycoerythrin)
Host species	Mouse
Conjugation	Phycoerythrin. Ex: 488nm, Em: 575nm
Tested applications	Suitable for: Flow Cyt
Species reactivity	Reacts with: Human
Immunogen	U937 cell line (Human).

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.1% Sodium Azide Constituents: 5% BSA, PBS
Purity	Ion Exchange Chromatography
Clonality	Monoclonal
Clone number	B-T47
Myeloma	x63-Ag8.653
Isotype	IgG1

Applications

Our [Abpromise guarantee](#) covers the use of **ab47210** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
-------------	-----------	-------

Flow Cyt

Application notes Flow Cyt: Use 10µl for 10⁶ cells or 100 µl of whole blood.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

Target

Function	Presents carbohydrate ligands to selectins. Also implicated in tumor cell metastasis.
Sequence similarities	Belongs to the LAMP family.
Post-translational modifications	O- and N-glycosylated; some of the 18 N-linked glycans are polylactosaminoglycans.
Cellular localization	Cell membrane. Endosome membrane. Lysosome membrane. This protein shuttles between lysosomes, endosomes, and the plasma membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors