

Product datasheet

# Anti-Streptococcal Pyrogenic Exotoxin C antibody ab53404

Overview

<b>Product name</b>	Anti-Streptococcal Pyrogenic Exotoxin C antibody
<b>Description</b>	Rabbit polyclonal to Streptococcal Pyrogenic Exotoxin C
<b>Host species</b>	Rabbit
<b>Specificity</b>	Minimal cross-reactivity with staphylococcal enterotoxins A through E, ET, TSST and alpha hemolysin and streptococcal Pyrogenic exotoxins A and B.
<b>Tested applications</b>	<b>Suitable for:</b> Indirect ELISA
<b>Species reactivity</b>	<b>Reacts with:</b> Streptococcus pyogenes
<b>Immunogen</b>	Streptococcal Pyrogenic Exotoxin C (SPEC)

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Preservative: None Constituents: 0.9% Sodium chloride, 10mM PBS
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab53404** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Indirect ELISA		Use at an assay dependent dilution.

## Target

---

### Relevance

Streptococcal Pyrogenic Exotoxin C is a causative agent of the symptoms associated with scarlet fever, has been associated with streptococcal toxic shock-like disease and may play a role in the early events of rheumatic fever.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors