

Product datasheet

Anti-ZNF259 antibody ab56419

2 Images

Overview

Product name	Anti-ZNF259 antibody
Description	Mouse monoclonal to ZNF259
Tested applications	Suitable for: WB, IHC-P
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment: IRELVTKNPF TLGDSSNPGQ TERLQEFSSQK MDQIIEGNMK AHFIMDDPAG NSYLQNVYAP EDDPEMKVER YKRTFDQNEE LGLNDMKTEG YEAGLAPQR, corresponding to amino acids 361-460 of Human ZNF259 Run BLAST with ExPASy Run BLAST with NCBI

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None PBS, pH 7.2
Purity	Protein G purified
Clonality	Monoclonal
Isotype	IgG2a
Light chain type	kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab56419** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		
IHC-P		

Application notes

WB: Use at a concentration of 1-5 µg/ml.

IHC-P: Use at a concentration of 3 µg/ml.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Target

Function

May be a signaling molecule that communicates mitogenic signals from the cytoplasm to the nucleus.


Sequence similarities

Belongs to the ZPR1 family.

Cellular localization


Cytoplasm. Nucleus > nucleolus. Translocates to the nucleolus after treatment with mitogens.

Images



IHC-P - ZNF259 antibody (ab56419)

ZNF259 antibody (ab56419) used in immunohistochemistry at 3µg/ml on formalin fixed and paraffin embedded human small Intestine.



Western blot - ZNF259 antibody (ab56419)

Predicted band size : 51 kDa

ZNF259 antibody (ab56419) at 1µg/lane + Jurkat cell lysate at 25µg/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors