

Product datasheet

Anti-C9orf61 antibody ab56533

1 Image

Overview

Product name	Anti-C9orf61 antibody
Description	Mouse monoclonal to C9orf61
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Rat, Human
Immunogen	Recombinant fragment: QIMAPLQPST SRAHRLPSRR QPGLLHLQSC GDLHTFTPAG RPRAERRPRR VEAERPHSLI GVIRETVL, corresponding to amino acids 383-451 of Human C9orf61 Run BLAST with ExPASy Run BLAST with NCBI

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None PBS, pH 7.2
Purity	Protein G purified
Clonality	Monoclonal
Isotype	IgG2a
Light chain type	kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab56533** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
-------------	-----------	-------

WB

Application notes WB: Use at a concentration of 1-5 µg/ml.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Target


Relevance

The function of C9orf61 is unknown.

Cellular localization

Single-pass type I membrane protein

Images



Predicted band size : 32 kDa

C9orf61 antibody (ab56533) at 1ug/lane +

PC-12 cell lysate at 25ug/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors