

Product datasheet

Anti-ZNF277 antibody ab56845

1 Image

Overview

Product name	Anti-ZNF277 antibody
Description	Mouse monoclonal to ZNF277
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Recombinant full length protein, corresponding to amino acids 1-279 of Human ZNF277

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None PBS, pH 7.2
Purity	Protein G purified
Clonality	Monoclonal
Isotype	IgG1
Light chain type	kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab56845** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		


Application notes WB: Use at a concentration of 1-5 µg/ml.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

Target

Function	May be involved in transcriptional regulation.
Sequence similarities	Belongs to the ZNF277 family. Contains 2 C2H2-type zinc fingers.
Cellular localization	Nucleus.

Anti-ZNF277 antibody images



Predicted band size : 52 kDa

ZNF277 antibody (ab56845) at 1ug/lane +
Jurkat cell lysate at 25ug/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors