

Product datasheet

Anti-ZNF277 antibody ab56845

1 Image

Overview

| | |
|----------------------------|---|
| Product name | Anti-ZNF277 antibody |
| Description | Mouse monoclonal to ZNF277 |
| Host species | Mouse |
| Tested applications | Suitable for: WB |
| Species reactivity | Reacts with: Human |
| Immunogen | Recombinant full length protein, corresponding to amino acids 1-279 of Human ZNF277 |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Storage buffer | Preservative: None PBS, pH 7.2 |
| Purity | Protein G purified |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Light chain type | kappa |

Applications

Our [Abpromise guarantee](#) covers the use of **ab56845** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-------------|-----------|---|
| WB | | Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 52 kDa. |

Target

| | |
|------------------------------|---|
| Function | May be involved in transcriptional regulation. |
| Sequence similarities | Belongs to the ZNF277 family. Contains 2 C2H2-type zinc fingers. |
| Cellular localization | Nucleus. |

Images



ZNF277 antibody (ab56845) at 1 ug/lane +
Jurkat cell lysate at 25ug/lane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors