

Product datasheet

Anti-UBL5 antibody ab57516

1 Image

Overview

<b>Product name</b>	Anti-UBL5 antibody
<b>Description</b>	Mouse monoclonal to UBL5
<b>Tested applications</b>	<b>Suitable for:</b> WB
<b>Species reactivity</b>	<b>Reacts with:</b> Human
<b>Immunogen</b>	Recombinant fragment: MIEVVCNDRL GKKVRVKCNT DDTIGDLKKL IAAQTGTRWN KVLKKWYTI FKDHVSLGDY EIHDGMNLEL YYQ, corresponding to amino acids 1-74 of Human UBL5 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Storage buffer</b>	Preservative: None PBS, pH 7.2
<b>Purity</b>	Protein G purified
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Light chain type</b>	kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab57516** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
-------------	-----------	-------

WB

**Application notes** WB: Use at a concentration of 1-5 µg/ml.

This antibody has only been tested in WB against the recombinant fragment used as

immunogen. We have no data on the detection of endogenous protein.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## Target


---

<b>Tissue specificity</b>	Ubiquitous. Highest level of expression in heart, skeletal muscle, kidney, liver, iris and lymphoblasts.
<b>Sequence similarities</b>	Contains 1 ubiquitin-like domain.
<b>Cellular localization</b>	Cytoplasm.

---

## Anti-UBL5 antibody images

---



Western blot against tagged recombinant protein immunogen using ab57516 UBL5 antibody at 1ug/ml. Predicted band size of immunogen is 34 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors