

## Product datasheet

# Anti-NTAN1 antibody ab67472

[1 Image](#)

### Overview

<b>Product name</b>	Anti-NTAN1 antibody
<b>Description</b>	Mouse polyclonal to NTAN1
<b>Tested applications</b>	<b>Suitable for:</b> WB
<b>Species reactivity</b>	<b>Reacts with:</b> Human
<b>Immunogen</b>	Full length NTAN1 protein (Human)
<b>Positive control</b>	NTAN1 transfected 293T cell lysate

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Storage buffer</b>	Preservative: None Constituents: Whole serum
<b>Purity</b>	Whole antiserum
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

Our [Abpromise guarantee](#) covers the use of **ab67472** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/1000. Detects a band of approximately 30 kDa (predicted molecular weight: 35 kDa).

### Target


<b>Relevance</b>	NTAN1 is required for the ubiquitin-dependent turnover of intracellular proteins that initiate with Met-Asn. This enzyme does not act on substrates with internal or C-terminal asparagines and
------------------	---

does not act on glutamine residues in any position.

## Cellular localization

Cytoplasmic

## Anti-NTAN1 antibody images



**All lanes :** Anti-NTAN1 antibody (ab67472) at 1/500 dilution

**Lane 1 :** NTAN1 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg/ml per lane.

### Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

**Predicted band size :** 35 kDa

**Observed band size :** 30 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors