

Product datasheet

Anti-ZNF557 antibody ab68259

1 Image

Overview

Product name	Anti-ZNF557 antibody
Description	Mouse polyclonal to ZNF557
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Full length ZNF557 protein (Human)
Positive control	ZNF557 transfected 293T cell lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: Whole serum
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab68259** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/1000. Detects a band of approximately 35 kDa (predicted molecular weight: 49 kDa).

Target

Function	May be involved in transcriptional regulation.
-----------------	--


Sequence similarities

Belongs to the krueppel C2H2-type zinc-finger protein family.
Contains 10 C2H2-type zinc fingers.
Contains 1 KRAB domain.

Cellular localization

Nucleus.

Images



Western blot - ZNF557 antibody (ab68259)

All lanes : Anti-ZNF557 antibody (ab68259)
at 1/500 dilution

Lane 1 : ZNF557 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at
1/2500 dilution

Predicted band size : 49 kDa

Observed band size : 35 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors